
eAdjust

ADAPTER

A.

B.

C.

D.

E.

F.

1 Information

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1 Information

Date of Last Update: April 30, 2024

- Please read this document carefully before using the product and observe the safety instructions.
- Contact the manufacturer if you have any questions about the product or if any issues arise.
- Report any serious incidents caused by the product, especially in cases of deteriorated health, to the manufacturer and the relevant authorities in the country concerned.
- Keep this document for future reference.

The product "eAdjust" is henceforth referred to as the product/ adapter.

This user manual provides you with important information on the use, adjustment, and handling of the product.

Only use the product in accordance with the information provided in the accompanying document.

2 Product Description

2.1 Function

The product is designed to be used as an alignment adapter with the capability to angle, offset, and rotate a custom-fitted transtibial (TT) or transfemoral (TF) prosthetic socket and is left in place in the finished prosthesis.

2.2 Compatibility

The product is designed to be used in combination with a standard 4-hole adapter manufactured for either laminated or thermoplastic of TT and TF prostheses. It can also be used in conjunction with the Direct Socket - standard 4-hole adapter. The pyramid on the eAdjust is of a standard type and fits with other manufacturers' female pyramid receivers.

2.3 Features

The adapter is made of grade 5 Titanium and 7075-T73 Aluminum. It has been tested according to ISO 10328 P6. This adapter allows for 23 mm horizontal slide, 11 degrees of angulation..

2.4 Lifespan

The adapter has been tested in accordance with ISO 10328 P6 for three million load cycles. This corresponds to a usage period of 3–5 years, depending on the user's level of activity.

3 Purposeful Use

3.1 Intended Use

The product is designed for use in the manufacture of prostheses for lower extremities. The adapter connects the prosthetic socket with other components in the leg prosthesis for amputation levels TT and TF, where the user has a maximum body weight of 125kg. The product has been developed and tested for everyday activities and must not be used for extraordinary activities or environments. Extraordinary activities include, for example, extreme sports (rock climbing, skydiving, paragliding, etc.).

WARNING: If the product is subjected to abnormal stress, accidents, or falls, the patient should contact their orthopedic engineer to have the product and its components inspected for possible damage

WARNING: If there is a change in, or loss of functionality, unusual noise, or if the product shows signs of damage or wear that impedes the normal function of the prosthesis, the patient should stop using the prosthesis and contact their Certified Prosthetist (CP).

3.2 Indications

- For users with transtibial (TT) or transfemoral (TF) amputation.
- For unilateral or bilateral amputation.

3.3 Contraindications

- Any conditions that contradict or exceed the instructions provided in the "Safety" and "Intended Use" sections.

3.4 Qualifications

The treatment of a patient with this product must only be carried out by an orthopedic engineer, or by qualified personnel under the direct supervision of an orthopedic engineer.

4 Safety

4.1 General Safety Instructions

If the safety instructions are not followed

injury to persons/products may occur if the product is used in certain situations.

- Follow the safety instructions and take the precautions specified in this accompanying document.

Use in the Presence of Children

Risk of injury if small parts are swallowed.

- Never leave children unsupervised during use.

Personal Manipulation of the Product

Risk of failure due to breaking of load-bearing parts or malfunction of the product.

- No work should be performed on the product other than what is described in this manual.
- This product must not be adjusted or modified by the patient/user or any other unqualified person.

Incorrect Alignment or Assembly

Risk of damage to prosthetic components.

- For proper and safe assembly, follow the instructions in this document and use the specified torques for tightening screws.

If the Specified Torque Settings for Screws are not Followed

Risk of failure due to breaking of load-bearing parts as screws may loosen.

- Follow the specified torque settings to ensure safe use.
- The specified torque settings must not be exceeded.
- Ensure there is no oil, dirt, or grease on the screws before thread locking is applied and the screws are tightened.

5 Assembly

5.1 Tools Required for Assembly

- Torque wrench
- Torx-T30
- Hex key 2.5mm
- Hex key 4mm

5.2 Assembly of Screw Fixture and Sleeve Adapter

1. Install the threaded screw fixture (1) into the prosthetic socket using an appropriate M6 screw. Apply Loctite 241 to the screws. Tighten them to **15Nm**.
2. Loosen the lock screw in the sleeve adapter. (4)
3. Install the sleeve adapter (4) onto the screw fixture. (1) **NOTE:** Ensure that the sleeve adapter is level with the screw fixture. Tighten the lock screw (5) to **8Nm** (see image A).

5.3 Setting the eAdjust Before Testing the Prosthesis

1. Ensure that the three lock screws (10) are not installed in the locking ring. (6)
2. Loosen the lock screw (7) in the locking ring slightly.
3. You can now slide/move the pyramid (8) to the desired position in preparation for testing the prosthesis on the patient.
4. Tighten the lock screw in the locking ring (7). **NOTE:** Tighten this to **8Nm** before the patient walks with the prosthesis. (see image B)
5. Install the prosthetic socket with the adapter onto other prosthetic components.

5.4 Adjustment of eAdjust During Prosthesis Testing

If the prosthetic socket needs adjustment during testing to optimize the setting of the prosthesis, this can be done easily without the patient having to remove the prosthesis or any other components needing to be disassembled.

1. It is recommended to initially mark the position and orientation of the locking ring with a pen. (see image C)
2. Loosen the lock screw in the locking ring until it can be slid by hand.
3. The adapter and socket can now be moved to the desired position. (see image C)
4. Tighten the lock screw (7). **NOTE:** Tighten this to **8Nm** before continuing with the prosthesis testing.

5.5 Locking the eAdjust Before Finalizing the Prosthesis

1. Mark the position and direction of the adapter to maintain the setting.
2. Unscrew the lock screw of the sleeve adapter (5) and apply Loctite 241 to the threads. Tighten to **8Nm**.
3. Unscrew the locking ring's lock screw and apply Loctite 241 to

the threads. Tighten to **8Nm**.

4. Remove the M10 stop screw (9) from the top of the pyramid.
NOTE: This should not be replaced. (see image D)
5. Turn down the M10 center screw (2) 5 turns. **NOTE:** It should be screwed clockwise from the top of the pyramid. (see image E)
6. Apply Loctite 241 as far down as possible on the inner thread of the pyramid.
7. Tighten the M10 center screw (2) (**NOTE:** counterclockwise) to a torque of **15 Nm**. Use a wrench or fixed key to hold the pyramid while the center screw is tightened with a torque wrench. (see image F)
8. Apply Loctite 241 to Lock screws (10) and tighten these to **1.5Nm**.

6 Maintenance

For user safety and to maintain operational reliability, we recommend annual safety inspections.

7 Cleaning

The product is water-resistant and washable. If exposed to saltwater, sweat, urine, or similar, it should be cleaned.

- Remove impurities from the product with a damp cloth and soap.
- Rinse the product with clean freshwater.
- Dry the product with a soft cloth.

WARNING: The product must not be stored in a damp condition. Allow any remaining moisture to air dry.

8 Waste Disposal

This product must not be disposed of with unsorted household waste. Improper waste management that does not comply with the regulations of your country can have harmful effects on the environment and health. Follow the disposal and recycling instructions provided by the relevant authority in your country.

9 Legal Information

All legal conditions are subject to the legislation in the country where the product is used and may therefore vary.

9.1 Liability

The manufacturer is only liable for damages that occur when the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer disclaims all liability for damages resulting from failure to follow the instructions in this document.

9.2 Warranty

The manufacturer's warranty for the product is valid from the date of purchase. The warranty covers defects that can be attributed to design faults, manufacturing faults, or material faults reported within the specified warranty period.

Warranty period: 12 months

9.3 CE Marking

Embreis AB hereby certifies that the product complies with applicable European regulations concerning medical devices.

Class 1 according to EU Regulation 2017/745

10 Technical Data

Product number: E0500HP

Weight: 258g

Connection: Pyramid

Maximum body weight: 125 kg

System height: 45 mm

Installation height: 33 mm

Material: Titanium, duralumin

Range of rotation: 360°

Angling: 11°

Horizontal slide: 23 mm

Lifespan: 5 years

SCAN ME

Embreis AB
Tumstocksvägen 11B
S-187 66 Täby
Sweden
Tel. + 46 8 410 621 00

info@embreis.com

www.embreis.com

UK Responsible Person
Wellkang Ltd
16 Castle St. Dover
Kent CT16 1PW,
England, United Kingdom

<http://www.wellkang.ltd.uk/>

embreis

Embreis AB | Tumstocksvägen 11B | S-187 66 Täby Sweden

Tfn + 46 8 410 621 00 | info@embreis.com

www.embreis.com