eSling

AXELSLINGA
SKULDERFATLE
SKULDERSLYNGE
SHOULDER SLING

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EN - eSling can be used on either the right-hand or the left-hand side. To change sides, simply pull the middle strap out from the buckle and turn it so that the stitched Velcro side of the strap is on the outside. This allows you to adjust the size during application. The actual Velcro should then be turned inward so that it can be fastened.

If you have any questions with regard to this product, please contact the party prescribed or dispensed the product to you, or get in touch with Embreis AB.

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1. Open up both the clip on the shoulder strap and the strap for the wrist and underarm. Place the underarm in the sling with the elbow resting on the rear edge.

2. Position the strap with the buckle at wrist level and tighten the sling. Place the other strap around the elbow and tighten.

3. Run the shoulder strap over the healthy shoulder and close the buckle, and centre the shoulder pad. Adjust the length of the shoulder strap by pulling on the end of the strap. Cut the shoulder strap to the desired length and place the Velcro at the end directly over the buckle.

4. Adjust the length and angle of the thumb grip by moving it to a suitable position on the shoulder strap or sling. The thumb grip is designed to provide comfort and relief for the wrist but can be easily removed if so desired.

5. Fasten the waist strap at wrist level. The length of the strap can be adjusted by removing the Velcro attachment at the end, cutting to the desired length and replacing the Velcro band.

6. Adjust the width of the waistband at the rear after fastening at the front of the sling. The length is adjusted by releasing the Velcro at the front of the sling and cutting it to the desired length. The waistband may be removed entirely by pulling it out of the sling.

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