### HINGE RANGE OF MOTION ADJUSTMENT

**A.** Brace is shipped with 10° hinge extension steps installed with additional flexion and extension steps in box. Steps are labeled in degrees:  
Extension steps: 0°, 10°, 20°, 30°, & 40°  
Flexion steps: 45°, 60°, 75°, & 90°  
(Brace is shipped with no flexion steps installed)

To change hinge steps, first remove the conyle pads held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

Next, use a screwdriver to remove screws. Remove existing step, insert new step, replace screw. Be sure that screws are seated and snug. Both medial (inside) hinge and lateral (outside) hinge steps must be identical in degree.

**OSTEOARTHRITIS (OA) HINGE INSTRUCTIONS**  
OA symptoms can be relieved by applying a lateral force to the outside of the knee by turning the circular dial in the hinge clockwise. With the knee bent in either a seated or standing position, turn the dial to the 0 to 10 degree increment. The dial will display settings ranging from 0 to 10 in 10 degree increments; the higher the number the greater the lateral force being applied. To reduce the lateral force, turn the dial counter-clockwise to a lower setting.

- Make sure the large gel conyle pad is attached to the lateral (outside) conyle basket.

### CARE AND CLEANING

Wash the brace with soap and water. Do not use bleach, disinfectant, or any other harsh cleaning agents. Wash the brace with mild detergent, rinse, AIR DRY. Occasionally apply a dry lubricant to hinge especially after use in water.

The frame pads are held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

The knee hinge conyle pads are held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.


To remove the brace straps first remove the strap pads on the inside of each strap. Then loosen the hook/loop tabs at each strap end and feed the straps through the anchor points on the side of the brace until the strap is free. Be sure to replace the straps and strap pads according to their strap number/location, see illustration.

### FITTING TIPS

To achieve the best possible fit of the brace straps, be sure that the strap pads on the inside of each strap are centered on the strap and provide a comfortable contact point between the skin and the strap.

Strap length adjustment: If a strap is too long or too short (preventing the numbered strap tabs from properly adhering to the strap), remove the strap pad and unfasten the strap tab, then shorten the "tail end" of the strap and slide the strap to either lengthen or shorten. Relasten the strap ends and replace the strap pads.

If there is no diagnosed or apparent instability associated with the accompanying knee OA, simply the step application for FUSION OA braces 4 and 6 in front of the knee. If there is a diagnosed or apparent instability associated with the knee OA, simply the step application for FUSION OA braces 4 and 6 in front of the knee.

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