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
CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER. FOR SINGLE PATIENT USE ONLY. DO NOT SHARE. DO NOT USE ON ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

CAUTION: CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO OPERATION. THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER FIT, POSITIONING AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.


ROADRUNNER KNEE BRACES

APPLICATION INSTRUCTIONS FOR KNEE BRACES:

The wrap-on brace is applied by pulling the brace on with the ratchet-coupler centered over the knee. Secure the closure firmly.

Both the wrap-on and pull-on braces have tightening straps. The straps should be tightened firmly, starting with the two straps closest to the knee.

The hinge settings are labeled in degrees as follows:
Extension settings: 0°, 10°, 20°, 30°, 40°, 45°, 60°, 75°, 90° & 110°
Flexion settings: 0°, 10°, 20°, 30°, 40°, 45°, 60°, 75°, 90° & 110°

To change the stops, first, completely remove the screw as shown. Then remove the stop. Insert the new stop along with the holes in the hinge. Then reinsert the screw. Using a Phillips head screw driver, ensure that the screws are seated in the recess and completely tightened. Both medial (inside) and lateral (outside) hinge stops must be identical in degree.

NOTE: Due to soft tissue compression, we recommend the use of extension stops at least 10° greater than the desired limit (e.g., if you desire a 10° extension limit, use the 20° extension stop).

- TIPS ON USE AND CARE OF YOUR KNEE BRACE:**
- After exposure to salt water or dirt, thoroughly rinse the brace with fresh water and allow to air dry.
 - A dry lubricant, such as Teflon spray, may be used on the hinge occasionally. Do not use wet lubricants.
 - The hinge assemblies may be removed to allow easier washing of the sleeve and straps.
 - The side bars may be bent to accommodate individual legs. Bend in amounts above and below the hinge.
 - Periodically inspect the sleeve, straps and hinges for excessive wear. Be sure all screws are tight and stop in place.

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