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

CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER. FOR SINGLE PATIENT USE ONLY. DO NOT SHARE. DO NOT USE ON ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

CAUTION: CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO OPERATION. THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.

ROADRUNNER KNEE BRACES


APPLICATION INSTRUCTIONS FOR KNEE BRACES:

The wrap-on brace is applied with the slip-on and pull-on braces by pulling the brace on with the ratchet-coutlet centered over the knee. Secure the closure straps. The slip-on brace is applied by pulling the brace on with the ratchet-coutlet centered over the knee. Secure the closure straps. Both the wrap-on and pull-on braces have tightening straps. The straps should be tightened firmly, starting with the two straps closest to the knee.

ADJUSTMENT INSTRUCTIONS

Polycentric hinge with stops

The hinge settings are labeled in degrees as follows:
 Extension stops: 0°, 10°, 20°, 30°, 40°
 Flexion stops: 45°, 60°, 75° & 90°
 (No flexion stops are installed upon shipment)

To change the stops, first, completely remove the screw as shown. Then remove the stop. Insert the new stop along with the holes in the hinge. Then, reinsert the screw. Using a Phillips head screw driver, ensure that the screws are seated in the recess and completely tightened. Both medial (inside) and lateral (outside) hinge stops must be identical in degree.

NOTE: Due to soft tissue compression, we recommend the use of extension stops at least 10° greater than the desired limit (e.g., if you desire a 10° extension limit, use the 20° extension stop).

TIPS ON USE AND CARE OF YOUR KNEE BRACE:

- After exposure to salt water or dirt, thoroughly rinse the brace with fresh water and allow to air dry.
- A dry lubricant, such as Teflon spray, may be used on the hinge occasionally. Do not use wet lubricants.
- The hinge assemblies may be removed to allow easier washing of the sleeve and straps.
- The side bars may be bent to accommodate individual legs. Bend in amounts above and below the hinge.
- Periodically inspect the sleeve, straps and hinges for excessive wear. Be sure all screws are tight and stop in place.

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