## **Xtern line** Sizing chart XTERNCLASSIC\* XTERNSUMMIT\* XTERNFRONTIER\*







SHOE SIZE		XTERN KIT CLASSIC			XTERN KIT SUMMIT			XTERN KIT FRONTIER		
us		PRODUCT: 100100-SA-L	PRODUCT: 100100-SA-M	PRODUCT: 100100-SA-S	PRODUCT: 101100-SA-L	PRODUCT: 101100-SA-M	PRODUCT: 101100-SA-S	PRODUCT: 102100-SA-L	PRODUCT: 102100-SA-M	PRODUCT: 102100-SA-S
MEN	WOMEN	LARGE	MEDIUM	SMALL	LARGE	MEDIUM	SMALL	LARGE	MEDIUM	SMALL
16		•			•			•		
15.5		•			•			•		
15								•		
14.5		•			•			•		
14								•		
13.5		•			•			•		
13								•		
12.5		•			•			•		
12		•			•			•		
11.5		•			0	•		•		
11		•			0	•		•		
10.5	11.5	0	•			•		0	•	
10	11									
9.5	10.5		•			•			•	
9	10									
8.5	9.5		•			•			•	
8	9		0	•		0	•		0	•
7.5	8.5		0	•		0	•		0	•
7	8			•			•			•
6.5	7.5			•			•			•
6	7			•			•			•
5.5	6.5			•			•			•
5	6			•			•			•
4.5	5.5						•			
4	5			•			•			•
3.5	4.5			•			•			•
3	4			•			•			•
2.5	3.5						•			

O Upgrade to alternate size according to shoe or foot width



Turbomed Orthotics
Turbomed Foot Drop Club (Customer club)



Turbomed Orthotics Inc. 240 Sydney, Suite 101 St-Augustin-de-Desmaures QC G3A 0P2 Canada

## Which Xtern is right for you?



## See through design that is almost invisible when worn

L-1951

**② 25%** lighter\*\*

L-1951

- ▼ 18% more dorsi-flexion power\*\*
- **35%** more medio-lateral strength\*\*

## Designed for patients with reduced dexterity and hand mobility such as stroke patients

L-1932 (pending)

- Unique front support design
- Easy rear entry thanks to revolutionary frontal design
- Magnetized velcro strap is easily attached with one hand
- **25%** more medio-lateral stability\*\*



**HCPCS PDAC approved**